


Solution de travail · Baño Base · Bagno Base · Roztwór Roboczy

ReplenisherRegenerator
Entretien
Reforzador
Integratore
Regenerator
Пополнитель**Starter****CD-S**
Art. No.
103354**Water**Wasser
Eau
Agua
Acqua
Woda
Вода**Working Solution**Arbeitslösung
Solution de Travail
Baño Base
Bagno Base
Roztwór Roboczy
Приготовление


500 ml

+




40 ml

+



460 ml

=




1000 ml

OR · ODER · OU · O · ALBO · ЛИБО

WORKING SOLUTION · ARBEITSLÖSUNG

OPTION 2

Solution de travail · Baño Base · Bagno Base · Roztwór Roboczy

WaterWasser
Eau
Agua
Acqua
Woda
Вода**CD-MR**conc.
Konz.**Starter****CD-S**
Art. No.
103354**Working Solution**Arbeitslösung
Solution de Travail
Baño Base
Bagno Base
Roztwór Roboczy
Приготовление

910 ml

+




50 ml

+



40 ml

=



1000 ml

OPTION 1 = Ansatz von Arbeitslösung aus Regenerator.

Mixing of working solution from replenisher.

Préparation de la solution de travail à partir de l'entretien.

El baño base se prepara a partir del reforzador.

Preparazione della soluzione di bagno base a partire dal integratore.

Przygotowanie roztworu roboczego z regeneratora.

OPTION 2 = Ansatz von Arbeitslösung direkt aus Konzentrat.

Mixing of working solution directly from concentrate.

Préparation de la solution de travail à partir du concentré.

El baño base se prepara directamente del concentrado.


Preparazione della soluzione di bagno base direttamente dalle concentrate.

Przygotowanie roztworu roboczego bezpośrednio z koncentratu.

RA-4 CD-MR**SP 108****COLOUR DEVELOPER REPLENISHER · FARBENTWICKLER-REGENERATOR**

Révélateur Chromogène Entretien · Revelador Reforzador Color · Integratore Sviluppo ·


Regenerator Wywoływacza Barwnego · Цветной Проявитель Пополнитель

for · für · pour · para · per · на · на **4 x 10 l** Art. Nr. 104435**REPLENISHER · REGENERATOR · ENTRETIEN · REFORZADOR · INTEGRATORE · REGENERATOR · ПОПОЛНИТЕЛЬ****Water**Wasser
Eau
Agua
Acqua
Woda
Вода**CD-MR**conc.
Konz.**Replenisher**Regenerator
Entretien
Reforzador
Integratore
Regenerator
Пополнитель900 ml
9 l

+

100 ml
1 l

=

1000 ml
10 l

15-30°C

REPLENISHMENT RATE · REGENERIERATE

Taux d'entretien · Tasa de Refuerzo · Tasso d'integrazione · Dawka Regeneracyjna · Нормы пополнения

108 ml/m²**TEMPERATURE / TIME · TEMPERATUR / ZEIT**

Température / Temps · Temperatura / Tiempo · Temperatura / Tempo · Temperatura / Czas · Температура / Время

38° C ± 0,3°C / 45 s